



**PANTHERA**  
SLEEP

# PANTHERA X3™

## NEXT GENERATION DORSAL FIN SLEEP APPLIANCE 100% NYLON CAD/CAM

SMALL, LIGHTWEIGHT, CUSTOM-MADE

MAXIMIZED COMFORT &  
PATIENT COMPLIANCE

BRUXISM &  
BREAKAGE RESISTANT

### NYLON POLYAMIDE TYPE 12


- Biocompatible, USP Level VI
- Light, smooth, flexible
- Strong & highly durable

### NO PREMATURE POSTERIOR CONTACT

- Additional TMJ protection

### UNDERCUT RETENTION

- No ball clasps



### EASY AND PATENT PENDING CLIP-ON TITRATION

- No tools required
- Precise & Secure
- Reduced risk of error
- Can move backwards 1mm

### TRIANGLE ANGLED AT 20°

- Allows drinking and talking


### ADJUSTABLE VERTICAL DIMENSION

- Without compromising oral cavity

US FDA 510(k) cleared | ISO 13485 : 2016 | CE | MDSAP



**EASY CARE**  
No storage in water




**COMPATIBLE**  
with the latest  
intra-oral scanning  
technologies



**3 YEAR**  
Warranty

# TITRATION

Simple and easy clip system




Backward movements  
are possible

Suggested starting position: 60%  
Possible to move back to 50%

## ANTERIOR OPENING AVAILABLE



## UPPER AND LOWER BANDS



## PLATEAUS

Tailored contact surfaces



CONTACT US: + 1 855 233-0388



WWW.PANTHERASLEEP.COM